



Release Notes TOM[®] Core Tray Organizing Manager Core

Web application: V. 3.0.15 (Sprint 46)

iOS application: V 1.2.2

Date: 02.04.2020



Table of Contents

- 1 New Features – iOS Application 3
 - 1.1 New handling of additional info 3
 - 1.2 Improved system feedback to the user 4
 - 1.3 Improved design for repair requests 5
 - 1.4 Link to validated methods for Aesculap products 6
- 2 New Features – Web Application 7
 - 2.1 Additional information in the Dashboard 7
 - 2.2 New specialties 8
 - 2.3 New Languages 8
- 3 Bug fixes – iOS Application 9
- 4 Bug fixes – Web Application 9
 - 4.1 Specialties order issue 9
 - 4.2 Wrong visualization in Dashboard for Items 9




1 New Features – iOS Application

1.1 New handling of additional info

In one of our last versions we launched the functionality to add information to the latest production. This was possible via opening the latest production and adding the information (including up to 3 pictures) there. Due to the fact that this was not intuitive and user friendly we changed the behavior.

- It is possible now to enter this information directly on the start screen or in the Assembling History.
- For all latest productions the icon ⓘ is shown. Select this icon and a dialog is opened which can be used to enter the additional information including the pictures directly. This saves time and only one step is needed to start the process.



- After the additional information has been added the icon changes to ⓘ indicating that for the marked production an additional information has been added. With selecting this icon the additional information could be reviewed and for the latest productions also changed. For older productions the dialog will be opened read-only.

TOM® Release Notes

1.2 Improved system feedback to the user

Since the last version it is only possible to open the latest productions. The system did not allow to open an older production but there wasn't any information shown to the user, why the production could not be opened anymore. To give feedback to the customer following has been implemented.

- The user could recognize which production is the latest one and which productions are old. This is visible due to a different style. The latest production is shown with a shadow and the older one without:




Figure 1: Latest production





Figure 2: Older production

- As it is visible above the info icon is missing for the older production. This shows that NO additional information has been added to the production and due to the fact that it is not the latest one the user could not add an information anymore.
- If the user tries to open an older production a dialog is shown to inform the user why it is not possible yet to open an older production. It is possible to open the additional information directly via this dialog.



TOM® Release Notes

1.3 [Improved design for repair requests](#)

The design of the repair requests has been improved. The UI is more intuitive and the information is presented to the user in a better way now:

14:23 Tue 31. Mar
Not Charging

Cancel
Create Repair Request

Description	Item code	Manufacturer	Missing in set	Set number	Speciality
DUROGRIP NADELHALTER HEGAR- MAYO 150MM	BM065R	Aesculap AG	Universal Set	XLS	UNKN

Select repair reason

- Surface discoloration
- Non-functional
- Deformation/Crack
- Corrosion
- REPLACED

Broken instrument at the handle..

33 of 250 characters used


+


➤ CREATE



1.4 [Link to validated methods for Aesculap products](#)

Via the website <https://instacount.aesculap.de/> it is possible to find all relevant documents and validated methods for Aesculap products. The link to this website is integrated into TOM now and the website will be opened showing the information for the selected item with the defined language directly.

- If the user opens an Aesculap item to check the details at the top of the screen the icon  is shown. If the user selects this icon the browser is opened and the mentioned website is shown:



- To navigate back to TOM the user has to use the back button of the iPad and to open TOM again. We will work on that to improve also this process.

2 New Features – Web Application


2.1 Additional information in the Dashboard

It was not visible in the Dashboard for the CSSD manager if a production has been packed completely or with missing items until now. Due to the fact that this is a very important information to check the quality of the CSSD process this is shown in the assembling history now.

- Three new columns are shown in the list of productions provide important information to the user. The CSSD manager is able now to recognize directly which of the productions have been packed incompletely and how many items were missing.

Set title >	Number >	Speciality	Requested / Actual Quantity	Missing Items >	Date v	Full name >	More info
Dental Set	1	DENT	9 / 9	–	02.04.2020, 11:38	Oliver Baumann	!
Micro Cervical Set	1	SPINE	22 / 21 ▲	-1	02.04.2020, 09:51	Oliver Baumann	–
Laparoscopic Set	1	LAPA	21 / 21	–	02.04.2020, 09:50	Oliver Baumann	–
Dental Set	4	DENT	28 / 25 ▲	-3	02.04.2020, 09:47	Oliver Baumann	–
Dental Set	1	DENT	9 / 9	–	02.04.2020, 09:11	Oliver Baumann	–
Test 02.04.2020	2	ANES / CALA	16 / 15 ▲	-1	02.04.2020, 09:11	Oliver Baumann	!

- Furthermore the CSSD manager could open all the additional information added to a production via the iPad packer app. This is possible with clicking on the ! in column “More Info”.



- All information will also be exported to the Excel file and the CSSD manager could use this for further analyses.


Note: In one of the next versions we will integrate further analyses about the status of productions.

TOM® Release Notes

2.2 [New specialties](#)

After feedback from customers we added two new specialties which could be chosen for a set. These are:

- Anesthesia ANES
- Cath-Lab CALA



2.3 [New Languages](#)

TOM is available in following new languages now:

- Vietnamese
- Polish

TOM® Release Notes

3 [Bug fixes – iOS Application](#)

- None

4 [Bug fixes – Web Application](#)


4.1 [Specialties order issue](#)

There was an issue for sets which were assigned to more than 1 specialty. It happened that the order of the specialties was reversed and this led to issues with the analyses.

- This has been fixed and the order of the specialties is always the same now.

4.2 [Wrong visualization in Dashboard for Items](#)

The visualization of the Items section in the Dashboard was wrong. The default pictures of the images were not shown correctly and there was a second scroll bar which did not work:



- This has been fixed and both issues have been removed:

