



Prontosan®
Product Range

NOTHING'S GONNA STOP YOU NOW.
INNOVATING WOUND BED PREPARATION

Prontosan® – Only a clean wound can heal

Wound bed preparation and infection prevention are a prerequisite

INDICATIONS

Prontosan® Wound Irrigation Solution and Prontosan® Wound Gel/Gel X are indicated for cleansing and moistening of acute, chronic, infected skin wounds, 1st and 2nd degree burns (also 3rd degree for Prontosan Wound Gel X). They prevent the biofilm formation.

Prontosan® Wound Irrigation Solution is also ideal for moistening encrusted dressings, or bandages prior to removal and for instillation in combination with negative pressure wound therapy. ⁽¹²⁾

Prontosan® Gel and Gel X act as an effective barrier to reduce microbial penetration through the dressing and to decontaminate the wound bed between dressing changes.


ADVANTAGES

- Management and prevention of biofilm reformation ^{(1) (2)}
- Helps to prevent infections ⁽³⁾
- Improved patient outcomes, including time to heal ⁽⁴⁾
- Well-known substances with low allergenic potential ⁽⁵⁾
- Can be used up to 8 weeks after first opening




THE OPTIMAL SOLUTION FOR THE MANAGEMENT OF BIOFILM

Conventional Wound Irrigation



Traditional applied irrigation solutions like 0.9% NaCl, water or Ringer's solution glide over the biofilm without removing it.


Prontosan®



The unique combination of Betaine Et Polyhexanide. Prontosan® physically removes debris, slough and bioburden.

In comparison with traditional wound irrigation solutions, wound treated with Prontosan® healed within a significantly shorter period of time. ^{(4) (6)}

PRONTOSAN® WAS SHOWN TO HAVE THE HIGHEST EFFICACY (IN VITRO TESTS) ⁽⁷⁾



Comparison of the efficacy of different wound irrigation solutions against biofilm (*Pseudomonas aeruginosa*) after 24 hours exposure time.

Prontosan® Wound Irrigation Solution achieved a significant reduction ($p < 0.001$) of the biofilm by 87%.

Modern wound bed preparation with soft debridement

Prontosan® Debridement Pad



NEW FAMILY MEMBER

The latest addition to the Prontosan® product family is the Prontosan® Debridement Pad. It has been designed to support the Wound Bed Preparation by providing rapid, effective and safe debridement of the wound. The pad frees the wound from coatings and debris and absorbs excess debris and slough thanks to its advanced microfibre technology.



ADVANTAGES

- Microfiber technology
- Rapid and soft debridement
- To obtain a clean wound bed and a clean skin surrounding
- Frees the wound from coatings and dead cell residues
- Absorbs excess exudates
- Intact tissue is spared

APPLICATION ROUTINE

Chronic wounds, granulating ones as well as critically colonized or infected wounds:



1 WOUND DEBRIDEMENT

Moisten Prontosan® Debridement Pad with Prontosan® Wound irrigation solution. Applying gentle pressure, wipe the wound in a circular motion.



2 RINSING


Following debridement, it is advisable to irrigate the wound thoroughly for a second time with Prontosan® Wound irrigation solution and let it soak if indicated.



3 LONG LASTING EFFECT


Prontosan® Gel and Gel X act as an effective barrier to reduce microbial penetration through the dressing and to decontaminate the wound bed between dressing changes.

Product Details

	Product	Pack Size	Product Code
	Prontosan® Solution		
	40 ml ampoule	24	400484
	350 ml bottle	10	400403
	1000 ml bottle	10	400446

Composition

0.1 % Polyaminopropyl Biguanide (Polihexanide), Betaine Surfactant

	Product	Pack Size	Product Code
	Prontosan Wound Gel	1	400515
	Prontosan Wound Gel X	1	400517
	Prontosan Wound Gel X	1	400508

Composition

0.1 % Polyaminopropyl Biguanide (Polihexanide), Betaine Surfactant, Purified Water, Glycerol, Hydroxyethylcellulose

	Product	Pack Size	Product Code
	Prontosan® Debridement Pad		
		3 Pads	3908456
		10 Pads	3908457

Composition

Cleansing sheet composed of polyester and polypropylene microfibre supported by a backing layer made from polypropylene.

This international brochure contains information which is targeted to a wide range of audiences and could contain product details or information otherwise not accessible or valid in your country.

B. Braun Medical AG | info.bbmch@bbraun.com | www.bbraun.com/prontosan

Lit. References:

- Efficacy of various wound irrigation solutions against biofilms. Seipp HM, Hofmann S, Hack A, Skowronsky A, Hauri A., ZfW 2005;4(5):160-163.
- Davis SC, Harding A, Gil J, Parajon F, Valdes J, Solis M Et Higa A "Effectiveness of a polyhexanide irrigation solution on methicillin-resistant Staphylococcus aureus biofilms in a porcine wound model" in International Wound Journal ISSN 1742-4801, 2017, 1-8, © 2017 Medicalhelplines.com Inc and John Wiley Et Sons Ltd doi: 10.1111/iwj.12734.
- Moore, M 0.1% Polyhexanide-Betaine Solution as an Adjuvant in a Case-Series of Chronic Wounds, Surg Technology International, 2016.
- Bellingeri, A. et al. "Effect Of A Wound Cleansing Solution On Wound Bed Preparation And Inflammation In Chronic Wounds: A Single-Blind RCT". Journal of Wound Care 25.3 (2016): 160-168. Web.
- Evaluation of the efficacy and tolerability of a solution containing propyl betaine and polihexanide. Romanelli M, Dini V, Barbanera S, Bertone MS. Skin Pharmacol Physiol 2010;23(Suppl 1):41-44.
- Kaehn, Eberlein - In vitro test for comparing for comparing the efficacy of wound rinsing solutions, British Journal of Nursing, 2009;11(11):S4-S10
- Siepp HM, Hofmann F, Hack A, Skowronsky Hauri A. Efficacy of various wound irrigation solutions against biofilms. Zeitschrift für Wundheilung 2005;4(5):160-163
- Instruction for use: Askina® Calgitrol® Ag, Askina® Calgitrol® THIN, Askina® Calgitrol® Paste.
- Opasanon S, Magnette A, Meuleniere F, Harding K. Askina® Calgitrol® Made Easy. Wounds International 2012; 3(1). Available from www.woundsinternational.com
- Trial C, Darbas H, Lavigne J-P, Sotto A, Simoneau G, Tillet Y, et al. Assessment of the antimicrobial effectiveness of a new silver alginate wound dressing: a RCT. J Wound Care. 2010 Jan;19(1):20-6.
- Wounds International. Using Askina® Calgitrol® Paste for the treatment of diabetic foot infection: case studies. London: Wounds International 2013. Available from www.woundsinternational.com
- Kim et al - Negative-Pressure Wound Therapy with Instillation: International Consensus Guidelines, PRSJournal 2013, DOI: 10.1097/PRS.0b013e3182a80586