Askina® Calgitrol® Ag

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INSTRUCTIONS FOR USE

1. PRODUCT DESCRIPTION
Askina® Calgitrol® Ag is a sterile dressing consisting of two layers:

- an absorbent polyurethane foam layer which provides for the absorption of exudate
- an ionic silver alginate matrix which provides for broad antimicrobial effectiveness and helps prevent contamination from external bacteria; in the presence of wound exudate, the silver alginate matrix helps maintain a moist wound environment conducive to natural healing conditions.

2. INDICATIONS
Askina® Calgitrol® Ag is indicated for the management of excised partial to full thickness wounds, stage I-IV pressure sores, venous ulcers, second degree burns and donor sites. It is also used for external use only and may be used in the management of infected wounds under medical supervision at the discretion of the physician.

3. PRECAUTIONS
In the event of clinical infection, topical silver does not replace the need for systemic therapy or other systemic antibiotics.

4. CONTRAINDICATIONS
Askina® Calgitrol® Ag should not be used on:

- Patients with a known sensitivity to silver
- Areas where the use of Askina® Calgitrol® Ag should not be used:
- Where the presence of the metal is contraindicated.

5. INFORMATION FOR USE
Askina® Calgitrol® Ag is very simple to apply. It is a sterile, pre-cut dressing which is equipped. Askina® Calgitrol® Ag may be left in place for up to seven days or until the dressing is saturated with wound exudate.

5.1 Skin Preparation
a) Cleanse the wound with sterile saline (0.9%) or an appropriate wound cleansing solution such as Protoson® antiseptic solution.

5.2 Dressing Application
Select the appropriate Askina® Calgitrol® Ag size that will completely cover the wound surface, ensuring a 3cm margin beyond the edges of the wound. If necessary, several dressings can be overlapped to cover very large wounds. Place the dressing with the silver matrix (dark surface) facing the wound. Secure Askina® Calgitrol® Ag with a secondary dressing or secure with surgical tape or bandage as appropriate.

5.3 Dressing Removal / Changes
Remove retention dressing if applicable. Gently lift Askina® Calgitrol® Ag from the wound. If the dressing is saturated, more dressing may be applied. Rinse with saline or Ringer solution. Follow procedures 5.1. and 5.2. to apply a new dressing.

5.4 SPECIAL NOTES
The dressing may initially appear to increase in size in the early stages of Askina® Calgitrol® Ag treatment. This is normal and occurs as any wound debris is removed from the edges of the wound. This clears the way for healing. In the management of moderately exuding wounds, Askina® Calgitrol® Ag can only make a limited contribution to the more conductive to healing. There are cases where healing is impaired as a result of underlying conditions. In these instances, Askina® Calgitrol® Ag alone will not be sufficient. Appropriate and suitable treatment of the underlying conditions will be necessary as well. Therefore, after 4 weeks of Askina® Calgitrol® Ag treatment, there has been no improvement then, in line with accepted wound management practice, the original diagnosis and overall therapy should be reassessed. Askina® Calgitrol® Ag should be left in place as long as possible in order to prevent trauma to the fragile newly formed tissue and to reduce cross contamination through frequent dressing changes. Thick necroses should be removed before applying Askina® Calgitrol® Ag. Clinicians' Healthcare Professionals should be aware that there is very limited data on prolonged use of the product. The product contains products and particularly in children and neonates. At present there are no contraindications for the pediatric use of silver. To avoid inappropriate treatment, it is recommended that official best practice guidelines for wound management be followed.

5.5 Russian Safety Standards (GOST-R)

7. PRESENTATION
Askina® Calgitrol® Ag is available in the following sizes:

- 10 x 10 cm • 15 x 15 cm • 20 x 20 cm

Warning: Re-use of single-use devices creates a potential risk to the patient or user. It may lead to contamination and/or impairment of functional capability. Contamination of the device may lead to injury, illness or death of the patient.

Sterile unless pack damaged or opened.

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