



ASKINA® SILNET

SOFT SILICONE WOUND CONTACT LAYER
B. Braun Wound Care. Right. From the start.

Askina® SilNet

Soft silicone wound contact layer



Non-adherent: protects fragile granulation tissue

Atraumatic dressing change

Soft and conformable

Allows passage of exudate through the dressing

Askina® SilNet is made of a conformable non-woven material, coated on both sides with soft silicone.


It is used as a wound contact layer, it conforms to the wound surface and adheres safely to the surrounding skin.

The wound site is protected from mechanical disruption during dressing changes. The use of Askina® SilNet thus minimizes the trauma associated with dressing change.

TIP FOR APPLICATION

Askina® SilNet can be used in combination with secondary dressings

- ▶ Askina® Foam
- ▶ Askina® Heel
- ▶ Askina® Pad




Indications

- ▶ Skin tears
- ▶ Traumatic wounds
- ▶ 2nd degree burns
- ▶ Fixation of skin grafts
- ▶ Donor sites
- ▶ Lacerations

- ① Remove the transparent liner.
- ② Apply Askina® SilNet on the wound surface and gently remove the blue liner.
- ③ Apply suitable absorbing secondary dressing, like Askina® Foam.
- ④ Mode of action of Askina® SilNet: its porous structure allows vertical passage of exudate into the secondary absorbing dressing, with no risk of maceration.

Soft silicone adhesive

Askina® SilNet is made of soft silicone, which provides safe and secure fixation, without stripping of epidermal cells during removal, unlike traditional acrylic adhesive.



Askina® SilNet	Size	Pcs/Pack	Article number
	5 cm x 7.5 cm	10	5195710
	5 cm x 7.5 cm	3	5195703
	10 cm x 7.5 cm	10	5197510
	10 cm x 7.5 cm	3	5197503
	10 cm x 18 cm	10	5191810
	10 cm x 18 cm	3	5191803
	20 cm x 30 cm	5	5192305

Manufactured by: B. Braun Hospicare Lt. | Collooney | Co. Sligo | IRELAND
www.bbraun.com

This document, its contents, including institutional data, information, trademarks and logos mentioned herein are the exclusive property of B.Braun. Any representation and / or reproduction, total or partial, of this document and its contents without the express prior consent of B. Braun, is strictly prohibited and constitutes an infringement of the intellectual property rights of B. Braun.
Non binding document and photographs