


# B. Braun Omnifix®

Three-piece single-use syringe



OEM – The B2B division of B. Braun.

# B. Braun Omnifix®

## Three-piece single-use syringe

### Product performance

- Highly transparent barrel
- High contrast black graduation
- Permanent marking
- Ideal readability
- Safe plunger backstop
- Plunger stopper with double sealing ring for slow aspiration and injection, qualified according to ISO 10993 resp. DIN EN 30993
- Qualified according to ISO 7886-1

### Product material

- Polypropylene
- Latex-free
- PVC-free

### Product sizes

- Colored (non-sterile): 3 - 20 ml
- Transparent (sterile and non-sterile): 3 - 30 ml

### Packaging

- Colored (non-sterile): PE bag, 1000 pcs.
- Transparent (sterile): Box, 100 pcs.
- Transparent (non-sterile): PE bag
  - » 3 ml: 5300 pcs.
  - » 5 ml: 3600 pcs.
  - » 10 ml: 2000 pcs.
  - » 20 ml: 1000 pcs.
  - » 30 ml: 800 pcs.

## Colored plunger

*new!*

### Luer slip cone

Size	Graduation	Cone	Color of plunger	Article No. non-sterile
3 ml	0.1 ml	centric	yellow	9200210
			red	9200220
			green	9200230
			blue	9200250
5 ml	0.2 ml	centric	yellow	9200510
			red	9200520
			green	9200530
			blue	9200540
10 ml	0.5 ml	centric	yellow	9201010
			red	9201020
			green	9201030
			blue	9201040
20 ml	1 ml	centric	yellow	9202010
			red	9202020
			green	9202030
			blue	9202040

## Transparent plunger

### Luer slip cone

Size	Graduation	Cone	Article No.	
			sterile	non-sterile
3 ml	0.1 ml	centric	4616025V	8508713
5 ml	0.2 ml	eccentric	4616057V	8502340
10 ml	0.5 ml	eccentric	4616103V	8502358
20 ml	1 ml	eccentric	4616200V	8508747
30 ml	1 ml	eccentric	4616308F	-



Luer slip



Luer lock

### Luer lock cone

Size	Graduation	Cone	Article No.	
			sterile	non-sterile
3 ml	0.1 ml	centric	4617022V	8508462
5 ml	0.2 ml	centric	4617053V	8508607
10 ml	0.5 ml	centric	4617100V	8508429
20 ml	1 ml	centric	4617207V	8508470
30 ml	1 ml	centric	4617304F	4617320

*new!*



*Three-piece syringes with colored plunger!*

### Minimizing medication errors

Medication errors are a serious safety problem. Color-coded syringes can help to ensure diversification.

Our recommendation for use:

Color of plunger	Field of application
Yellow	Epidural anesthetics
Blue	Standard solutions
Red	Muscle relaxants
Green	Your drug portfolio

*Tailor-made solutions for the  
healthcare industry.*

*If you cannot find a suitable product in our  
product range, we are happy to provide you  
with your individual solutions – from product  
variation to complete new developments.*

*Your demand is our challenge:  
Whatever you need, we have the solution.*

*Headquarters:*

**B. Braun Melsungen AG**

OEM Division

Carl-Braun-Strasse 1

34212 Melsungen

Germany

Phone: +49 5661 71 4583

[www.bbraunoem.com](http://www.bbraunoem.com)

[info-europe@bbraunoem.com](mailto:info-europe@bbraunoem.com)

Edition: 01-2014