

SOLUTION FOR UROSTOMATES

YOUR NEEDS, OUR PRIORITIES

CREATE YOUR PATIENT SOLUTION

ALL STOMA AND DISTAL B. BRAUN BAGS ARE COMPATIBLE*

SKIN PROTECTOR ATTACHED
OR INDEPENDENT OF
THE COLLECTION BAG


SEVERAL FEATURES TO MAKE
B. BRAUN'S SOLUTION EASIER

COMPLETE RANGE OF ADULT
URINE COLLECTION BAGS
INCLUDING LEG AND BEDSIDE
BAGS

BAGS


1-PIECE SYSTEM

FLEXIMA®



URO SILK FLAT URO SILK CONVEX

FLEXIMA® ACTIVE



ACTIVE O'CONVEX URO

2-PIECE SYSTEM

FLEXIMA® 3S URO



FLAT CONVEX

PROXIMA® 2 URO



BORDER STANDARD CONVEX

FEATURES

UNIVERSAL CONNECTION

PROXIMA® 2 URO




ANTI-REFLUX VALVE



COLLECTOR BAGS

500 ML, 800 ML, 1,5L & 2L

B. BRAUN URIBAG



URIMED® B'BAG



URIMED® TRIBAG PLUS



*use our urostomy bag connector

STRENGTHS OF B. BRAUN SOLUTION

2

☀️ OUTLET FOR DAY

- Position open/closed
- Fits in perfectly strengthened by a click and a final locking system
- Easy handle; no risk of wet fingers

🌙 OUTLET FOR NIGHT

- Click at connection
- Can be connectd to a night drainage bag
- Outlet designed to avoid bag twist



CONVEXITY

- Low risk of leakage episodes
- Protection against peristomal skin disorder

COMFORTABLE FOR DAILY ACTIVITIES

- Universal connection
- Non kinking tube & large diameter

BENEFITS OF THE B. BRAUN SOLUTION

3

PATIENTS: PHYSICAL AND PSYCHOLOGICAL CARE

Security and easiness


No residual drop
Prevents urine backflow
Easy to handle

Comfort

Recovery during sleep
9 out of 10 urostomates have trouble sleeping¹

Autonomy

Physical, social and mental wellbeing



NURSES: ENSURE A GOOD EDUCATION OF PATIENTS

FACILITATES LEARNING OF THE DEVICE

- Dexterity disorder
- Anxiety, refusal of the disease

HELPS TO ACCEPT THE DISEASE

- No risk of leaks with the 2L collector bag
- Gives the possibility to make full nights

FACILITATES ORGANIZATION OF THE URINARY MONITORING IN THE HOSPITAL

- Approximative graduation, extension of the spare time interval
- Evaluation of urine quality and quantity

1. Cavdar I et al. Sleep and quality of life in people with ileal conduit. Scand J Urol 2016;50(6):472-476

BAGS

FLEXIMA® URO SILK FLAT – BOX OF 30*

Size	Ø(mm)	References ⁽¹⁾		Volume
Midi	✂12-55	● 044914	○ 044913	±630ml
	20		○ 44820	
	25		○ 44825	
	30		○ 44830	
	35		○ 44836	








*with 30 universal connectors

FLEXIMA® URO SILK CONVEX – BOX OF 10*

Size	Ø(mm)	References ⁽¹⁾		Volume
Midi	✂15-25		○ 44915	±615ml
	✂15-35	● 44919	○ 44917	
	25		○ 44920	
	30		○ 44922	

*with 10 universal connectors

+ BAGS

	References ⁽¹⁾	Quantity
	B. Braun Uribag - collection bag 2L 5071	30
	Urimed® Tribag Plus - 800L 3-Chamber leg bag 28306	10
	Urimed® Bag Plus - Leg bag 500mL 28501	10
	Urimed® B'Bag Sterile bags with NRV* - 2L - 90cm 29440	100
	Urimed® B'Bag Non sterile bags with NRV* - 2L - 90cm 29450	200
	Urimed® B'Bag Non sterile bags with NRV* - 1,5L - 90cm 29455	200
	Urimed® B'Bag Non sterile bags with NRV* - 2L - 90cm hanging accessory 29451	30 + 1

*NRV = Non-Return Valve

✂ Cut-to-fit | ● Beige | ○ Transparent | ● Split Cover / Beige with inspection window | *NRV = Non-Return Valve

(1) For most countries, a specific letter is added at the end of the radical reference

FLEXIMA® ACTIVE O'CONVEX URO – BOX OF 10*

Size	Ø(mm)	Quantity	References ⁽¹⁾	Volume
Midi	✂15-35	Box of 10	○ 4641715	±670ml

*with 10 universal connectors

FLEXIMA® 3S URO – BOX OF 30*

Size	Ø(mm)	References ⁽¹⁾		Volume
Midi	45	● 934045	○ 934145	±675ml
	55	● 934055	○ 934155	±635ml
	65	● 934065	○ 934165	±620ml






*with 30 universal connectors

PROXIMA® 2 UROSTOMY – BOX OF 10*

Size	Ø(mm)	References ⁽¹⁾	Volume
Midi	40	○ 73540A	±670ml
Midi	50	○ 73550A	±660ml
Midi	60	○ 73560A	±650ml

*with 5 universal connectors

+ ACCESSORIES

	References ⁽¹⁾	Quantity
	Superfiller Paste F05008	Tube 60 g
	Urimed® Fix - Leg fixation of urine bag 68520 (small)	2
	68530 (medium)	2
	68540 (large)	2
	Urimed® Straps - Leg fixation of urine bag 68550	2
	Urimed® Handle - Bag hanger at bedside or wheelchair 68548	1
	Urimed® Klett - Fixation of tubes and catheters 28260 (small)	10
	28261 (large)	10

This document, its contents, including institutional data, information, trademarks and logos mentioned herein are the exclusive property of B. Braun. Any representation and / or reproduction, total or partial, of this document and its contents without the express prior consent of B. Braun, is strictly prohibited and constitutes an infringement of the intellectual property rights of B. Braun.

Commercialized and manufactured by :

B. Braun Medical SAS | 26 rue Armengaud | 92210 Saint-Cloud - France

Tél. 01 41 10 53 00 | Fax 01 70 83 45 00 | www.bbraun.fr | www.bbraun.com

Société par Actions Simplifiée with a share capital of 31 000 000 € | RCS Nanterre 562050856

All stoma bags volumes are measured according to ISO 8670-2 standard.

This medical device is a regulated health care product and on the basis of this regulation is certified CE.

The products commercialized by B. BRAUN MEDICAL shall exclusively be used in accordance with the package inserts which have been updated and included in the boxes.

Non binding documents and photographs

For healthcare professional use only.