

Alcohol Pads B. Braun | Technical Data Sheet

Pads for Skin Cleansing – Isopropyl Alcohol 70% (non sterile)



	Alcohol Pads B. Braun		
REF	9160612		
PAD SIZE	Folded	32.5 mm x 30 mm	
	Unfolded	65 mm x 30 mm	
INGREDIENTS (INCI NAME)	Isopropyl Alcohol, Aqua		
MATERIALS/COMPOSITION	Solution of isopropanol	2-Propanol	70 percent by volume
		Purified water	ad. 100 volume percent
	Each pad is soaked with 0.62 g (Δ 0.7 ml) of an aqueous solution of isopropanol (70 vol. %) equivalent to 0.39 g (Δ 0.49 ml) isopropanol		
	Wipe	Polypropylene	100 percent by weight
		No chemical binders	
LOGISTICS	Packaging	Individually foil-wrapped, air-tight seal	
		Box with 100 pieces	
	Shelf life	4 years	
STERILIZATION	Non sterile		
APPLIED STANDARDS	Regulation EC 1223/2009 on cosmetic products		
	ISO 22716:2007 Cosmetic guidelines on good manufacturing practices (GMP)		
	ISO 9001		
PRODUCT CLASSIFICATION	Cosmetic product		

Abbreviation INCI International Nomenclature of Cosmetic Ingredients