

**NEW**



Aesculap Neurosurgery

# Aesculap® KOPITNIK AVM Micro Clip System

sterile packed | steril verpackt

**Aesculap provides all aneurysm clips sterile.**

From now on AVM Micro Clips are also available in a sterile packaging.

**Advantages of sterile packed AVM clips are...**

- ...ready to use
- ...no reprocessing costs

**Aesculap liefert alle Aneurysmen Clips steril.**

Von nun an sind auch die AVM Mikroclips in einer sterilen Verpackung verfügbar.

**Die Vorteile von steril verpackten AVM Clips sind...**

- ...sofort implantierbar
- ...keine Aufbereitungskosten

**Curved clips are only available at Aesculap!**  
**Die gebogenen Clips, exklusiv nur bei Aesculap!**

Art. No.		closing force between 50 g and 70 g Schließkraft zwischen 50 g und 70 g		Art. No.		closing force between 50 g and 70 g Schließkraft zwischen 50 g und 70 g	
Single packed		Blade length (mm)	Maximal opening (mm)	Single packed		Blade length (mm)	Maximal opening (mm)
Art. Nr. Einzelverp.		Mauillänge (mm)	maximale Öffnung (mm)	Art. Nr. Einzelverp.		Mauillänge (mm)	maximale Öffnung (mm)
FE952K		2	1.5	FE962K		2	1.5
FE953K		3	1.7	FE963K		3	1.7
FE954K		4	2.0	FE964K		4	2.0
FE955K		5	2.3	FE965K		5	2.3