


Omnican® Lance soft

Sterile single use lancets for capillary blood collection



Omnican® Lance soft				
Intended use		Lancets for capillary blood collection to be used with compatible lancing devices		
Compatible lancing devices		B. Braun Omnican® Lance, B. Braun Omnilance® and other lancing devices*		
REF		9167536E	9167510E	9167554A
Sales unit		200 pcs.	200 pcs.	50 pcs.
Languages		bg, cs, de, en, fr, it, nl, pl, ro, sk, sr, tr		en, id, ms, vi
Needle size	Gauge	28 G	30 G	30 G
	Diameter	0.36 mm	0.30 mm	0.30 mm
	Bevel	three-facet grinding		
	Puncture depth	variable acc. to lancing device used 7 different puncture depths with Omnican® Lance lancing device 5 different puncture depths with Omnilance® lancing device		
Logistics/ Packaging	Sales unit	200 pcs.	200 pcs.	50 pcs.
	Transport case	21 x 200 pcs.	21 x 200 pcs.	24 x 50 pcs.
	Shipping box	84 x 200 pcs.	84 x 200 pcs.	192 x 50 pcs.
Material	Needle	stainless steel, silicone oil coated		
	Lancet cap and body	Low Density Polyethylene (LDPE)		
	Latex-free	yes		
Shelf life		5 years		
Sterilization		Gamma ray irradiation		
Certificates		EN ISO 13485		
Product classification		Medical device acc. to MDD 93/42/EEC, CE class II a		



* e. g. from Abbott, Bayer, LifeScan, Menarini, Owen Mumford, Sanofi Aventis and Ypsomed.

All data refer to typical single values and may be subject to alterations.