

MANUFACTURER DECLARATION

Cleaning and disinfection of Omnitest® monitors in multi-patient use



In the context of multi-patient use of Omnitest® blood glucose monitors (one device is used for several patients), regular cleaning and disinfection of the device is recommended. Especially devices that have come in contact with blood inherit a potential risk to transmit infectious diseases. However, it is important to note that some cleaning and disinfection liquids can damage the outer surface of the Omnitest® devices.

WE HEREWITH CONFIRM THAT

Meliseptol® Wipes sensitive and 0.5 % ready-to-use
Hexaquart® plus solution are suitable for cleaning and
disinfection of Omnitest® devices.



Based on a very strong rationale the likelihood of a disease transmission by a "non-critical" medical device is far lower if the device is cleaned and disinfected using VAH-listed Meliseptol® Wipes sensitive or alternatively 0.5 % ready-to-use Hexaquart® plus solution compared to just cleaning.

REMARK:

Only damp wipe the outer surface of the device, do not allow liquids to enter the unit housing during cleaning and disinfection!

As designed for patient safety and environmental protection, please follow this recommendation carefully. In case of defected Omnitest® devices despite correct handling please get in contact with your local B. Braun representative for device exchange.

B. Braun Melsungen AG
CoE Diabetes Care

Felix Kever
International Product Manager

B. Braun Medical AG
CoE Infection Control

Alexander Müller
International Product Manager