

Urimed® Vision

Connect to freedom



Ordering information

Urimed® Vision Standard

Description	Ø 25	Ø 29	Ø 32	Ø 36	Ø 41
Article number*	IH2525	IH2529	IH2532	IH2536	IH2541
Units/box	30	30	30	30	30

Urimed® Vision Short

Description	Ø 25	Ø 29	Ø 32	Ø 36	Ø 41
Article number*	IH4525	IH4529	IH4532	IH4536	IH4541
Units/box	30	30	30	30	30

Urimed® Vision Ultra

Description	Ø 25	Ø 29	Ø 32	Ø 36	Ø 41
Article number*	IH3525	IH3529	IH3532	IH3536	IH3541
Units/box	30	30	30	30	30

Urine Bags

Description	Urimed Bag Plus 500ml	Urimed Tribag Plus 800ml	Urimed Bag Plus 1500ml	Urimed Bag 2000ml
Article number*	28501	28306	28150	28300
Units/box	30	10	30	30

*For each country a specific letter is added at the end of the reference radical.

Male external catheters

B. Braun Medical S.A.S. | RCS Nanterre 562050856 | 204, avenue du Maréchal Juin
92660 Boulogne Billancourt | France
www.bbraun.com

For healthcare professional use only - See instructions for use for contraindications and precautions.
Reference : ZG00713 - SOUS LE VENT - 05 61 63 93 68 - Fotolia

The content of this document, in particular the data, the information, the trademarks and the logos are B BRAUN MEDICAL SAS's sole property and are protected under copyright ©2010 BRAUN MEDICAL, all rights reserved. Consequently, all representation and/or reproduction, whether in part or in full, is forbidden and would be considered a violation of B BRAUN MEDICAL SAS's copyright and other proprietary rights. The products commercialized by B BRAUN MEDICAL shall exclusively be used in accordance with the package inserts which have been updated and included in the boxes.

B | BRAUN
SHARING EXPERTISE

Urimed® Vision: close to healthy skin






Male external catheters offer a discreet and reliable solution to male urinary incontinence. The real benefit for patients is to continue having an active life with minimal inconvenience in a secured way.


Urimed® vision is B. Braun male external catheter, a complete range made of silicone material that respects the skin properties and makes it comfortable to wear.

The **silicone** is a high performance, hypoallergenic and soft material.



-  **Moisture vapour transmission excellent gas exchange⁽¹⁾**
-  **High performance in water vapour permeability⁽¹⁾ compared to latex and elastomer**
-  **Transparent colour for a good monitoring of skin condition**

High performance in water vapour permeability⁽¹⁾



Urimed® vision offers to patients...

Security

Skin friendly

Kink proof silicone funnel helps to ensure uninterrupted urine flow
Secured connection to all major types of collecting bags.

Autonomy

Connected to urine bag, Urimed® leg & bed bags give patient **24 hours of autonomy** and freedom.

Ease of use

One piece system easy to apply and to remove, secured with an acrylic adhesive that fixes it to the skin.

A perfect fit

5 different diameters and 2 lengths.
Standard version for most patients and the **short one** for retracted penis or paediatric use.

High permeability for low risk of maceration



Choose the right connection !

MEC + Urine bag = **24h⁽²⁾ autonomy**

(1) Internal method according to the NF ISO 2528 Standard.
 (2) Refers to Instructions for use for a daily change.