




# Infusomat<sup>®</sup> plus Line Portfolio

## Product Information

Infusomat <sup>®</sup> plus Line portfolio	Description	Tube material	Length (cm)	Priming vol. (ml)	Y-port	Code no. (REF)
<b>Enteral Nutrition – ENFit™ Lines</b>						
	<ul style="list-style-type: none"> <li>Inline drip chamber</li> <li>Different tube materials</li> <li>Different lengths</li> <li>ENFit™ patient connector</li> </ul>	PUR**	330/220	28.0	needle-free	8700370
	<ul style="list-style-type: none"> <li>ENPlus Spike for enteral nutrition diets with ENPlus (cross) port</li> </ul>	PVC*	330/220	23.1	needle-free	8700400
	<ul style="list-style-type: none"> <li>Stopcock with male ENFit™ connector for enteral syringe with female ENFit™ connector</li> </ul>	PVC*	260/150	18.5	needle-free	8700380
<b>Regional Anesthesia – NRFit™ Line</b>						
	<ul style="list-style-type: none"> <li>15 µm fluid filter</li> <li>Drip chamber incl. AirStop function</li> <li>NRFit™ patient connector</li> <li>translucent yellow tube</li> </ul>	PUR**	300/200	21.8	–	8700410

\* PVC = DEHP-free; \*\* PUR = PVC-/DEHP-free



# Infusomat<sup>®</sup> plus Line Portfolio






## Product Information

# Infusomat<sup>®</sup> plus Line Portfolio





## Product Information

### Time to Care

The color codes and shapes of Infusomat<sup>®</sup> compact<sup>plus</sup> pump make it easy to insert the Infusomat<sup>®</sup> plus lines. The dedicated tubing material allows continued accuracy over long periods (96 h). The Anti Free-Flow Clamp automatically closes the line to prevent any risk of free-flow when the door is opened. In addition, the graphical user interface allows intuitive operation and guided loading reduces training needs.

Infusomat <sup>®</sup> plus Line portfolio	Description	Tube material	Length (cm)	Priming vol. (ml)	Y-port	Code no. (REF)
<b>SafeSet Lines</b>						
	▪ 15 µm fluid filter	PVC*	240/150	17.3	–	8700200
	▪ Drip chamber incl. AirStop function	PVC*	300/200	21.2	–	8700210
	▪ Different tube materials	PUR**	240/150	17.8	–	8700220
	▪ Different lengths	PVC*	300/200	22.2	needle-free	8700240
	▪ Luer-Lock Patient connector	PUR**	300/200	22.9	needle-free	8700250
<b>Standard Lines</b>						
	▪ 15 µm fluid filter	PVC*	240/150	17.3	–	8700310
	▪ Spike drip chamber	PUR**	240/150	17.8	–	8700320
	▪ Different tube materials	PVC*	300/200	22.2	needle-free	8700330
	▪ Different lengths	PVC*	300/200	22.2	needle-free	8700330
▪ Luer-Lock patient connector	PVC*	300/200	21.7	needle-based	8700340	
<b>Transfusion Lines</b>						
	▪ 200 µm blood filter	PVC*	240/150	17.3	–	8700350
	▪ Spike drip chamber					
	▪ Different lengths					
▪ Spinlock patient connector	PVC*	300/200	21.9	needle-free	8700360	
<b>UV-light protect Lines</b>						
	▪ 15 µm fluid filter	PUR**	240/150	17.8	–	8700260
	▪ Drip chamber incl. AirStop function					
	▪ Different lengths					
▪ Luer-Lock patient connector	PUR**	300/200	22.9	needle-free	8700270	
<b>Oncology Line</b>						
	▪ 15 µm fluid filter ▪ 0.2 µm hydrophilic fluid filter (Sterifix <sup>®</sup> ) ▪ Drip chamber incl. AirStop function ▪ Eurofix Spinlock patient connector	PUR**	240/150	20.6	–	8700300

\* PVC = DEHP-free; \*\* PUR = PVC-/DEHP-free

Infusomat <sup>®</sup> plus Line portfolio	Description	Tube material	Length (cm)	Priming vol. (ml)	Y-port	Code no. (REF)
<b>Cyto-Set<sup>®</sup> Lines</b>						
	<b>Standard Lines:</b>					
	▪ 15 µm fluid filter	PUR**	210/155	23.9	3 needle-free valves	8700480
	▪ Grip plate with integrated back-check valves					
	▪ Cytostatic resistant material	PUR**	210/155	23.9	5 needle-free valves	8700490
▪ Eurofix Spinlock patient connector						
▪ 0.2 µm hydrophilic fluid filter (Sterifix <sup>®</sup> )	PUR**	225/170	27.4	5 needle-free valves	8700470	
	<b>AirStop Lines:</b>					
	▪ Drip chamber incl. AirStop function	PUR**	210/155	23.9	3 needle-free valves	8700420
	▪ Grip plate with integrated back-check valves					
	▪ Cytostatic resistant material	PUR**	210/155	23.9	5 needle-free valves	8700430
▪ Eurofix Spinlock patient connector						
▪ 0.2 µm hydrophilic fluid filter (Sterifix <sup>®</sup> )	PUR**	225/170	27.4	5 needle-free valves	8700460	
	<b>UV-light protect Lines:</b>					
	▪ Light protection up to 520 nm	PUR**	210/155	23.9	3 needle-free valves	8700440
	▪ Drip chamber incl. AirStop function					
	▪ Grip plate with integrated back-check valves					
▪ Cytostatic resistant material	PUR**	210/155	23.9	5 needle-free valves	8700450	
▪ Eurofix Spinlock patient connector						
<b>Flushing Set</b>						
	▪ 15 µm fluid filter ▪ Inline drip chamber incl. AirStop function ▪ Spinlock patient connector	PUR**	300/200	22.4	needle-free	8700280

\* PVC = DEHP-free; \*\* PUR = PVC-/DEHP-free