

# Omnitest<sup>®</sup> 5

## MY SMART PARTNER

# Live your life with diabetes

Keep living your life the way you want to with Omnitest® 5 and the Om

- ✔ You can measure your values easily with Omnitest® 5 and obtain reliable data
- ✔ Your readings are automatically transferred to the Omnitest® app

HIGH ACCURACY  
ACCORDING TO  
EN ISO 15197: 2015

 **Bluetooth®**


### Easy data transfer

Your meter presents reliable readings and transfers them automatically to the Omnitest® app.



nitest® app as smart partners on your side.

- ✓ The digital diary helps you to analyze data and understand your diabetes
- ✓ Everything at hand without extra effort



### Approved quality

As one of the first diabetes apps to be awarded the test seal from DiaDigital, the Omnitest® app was praised for its ease of use and accessible design.

# Understand and manage your diabetes

## B. Braun's free Omnitest® app allows you to keep track of your treatment

- ✓ Enter your blood glucose readings manually or receive them automatically from your Omnitest® 5 meter with Bluetooth function
- ✓ Additional input options for diabetes-related data help you to better understand what influences your glucose levels

### THREE STEPS

#### Connect with Omnitest® 5

1



Go to the settings tab on the app and click on "Add new meter" in your app. Press the up and down buttons on your Omnitest® 5 meter simultaneously for at least 3 seconds.

2



Press "scan for meter" in your app to identify your meter. The app starts searching for Bluetooth devices.


3



Done! Once paired, all available data is immediately transferred to the Omnitest® app. From now on, new readings will be automatically transferred straight away.

### COMFORTABLE INPUT

#### Record your data



# tes more easily

-related data – anywhere, anytime.

- ✓ Clear graphical views and statistical modes simplify all given information and make it easy to analyze – you can recognize patterns and better understand your diabetes
- ✓ And if you want to share your information with others, you can easily export your data

## FULL ANALYSIS

Understand your diabetes



# Reliable readings you can count on

Use Omnican® Lance and Omnican® Lance soft lancets for smooth blood sampling.

## Virtually pain-free

- ✔ A gentle puncture and an ideal blood flow
- ✔ The lancet device allows for a straight and fast lancet leading

## Easy to operate

- ✔ The color signal upon loading shows when the device is ready
- ✔ One-hand operation to load and release




Image in original size

With its reliable measurement results Omnitest® 5 delivers the solid basis for your treatment.

## Comfortable

- ✓ Test strip port and screen are illuminated – take measurements whenever you need to
- ✓ Receive your measurement result and its deviation from your target level

## Preventive

- ✓ Automatic code recognition eliminates potential user errors
- ✓ Convenient test strip ejector allows for safe disposal of the used test strip



Large, illuminated test strip port

Convenient strip ejector

Large, illuminated screen

User-friendly buttons

Image in original size

# Omnitest® 5 with Bluetooth

## TECHNICAL INFORMATION

Sample volume	0.5 µL
Measuring time	5 seconds
Measuring range	10 – 600 mg/dL
Hematocrit range	20 – 60 %
Calibration	Plasma
Enzyme	FAD glucose dehydrogenase
Operating temperature	10 °C – 40 °C or 50 °F – 104 °F
Operating humidity	10 – 90 %
Screen type	Backlit LCD
Dimensions	78 mm x 43 mm x 16 mm (L x W x H)
Weight	45 g (including batteries)
Power supply	Two 3V lithium batteries (CR2032)
Data transfer via Bluetooth® technology	Omnitest® app
More info:	<a href="http://www.bbraun.com/diabetes-datamanagement">www.bbraun.com/diabetes-datamanagement</a>

## PRODUCT OVERVIEW

### Omnitest® 5

Set black	(mg/dL)
Set white	(mg/dL)
Solo white	(mg/dL)
Test strips	(2 x 25 pc or 2 x 50 pc)
Control solution	(M + H)
USB cable	(standard micro-USB)

### Omnican® Lance

Lancing device

### Omnican® Lance soft


Sterile lancets 28G or 30 G (200 pc)

### Omnitest® 5 Set contains:

blood glucose meter incl. 2 batteries, check strip, meter case, 10 test strips, lancing device, 10 lancets

### Omnitest® 5 Solo contains:

blood glucose meter incl. 2 batteries, check strip, meter case



B. Braun Melsungen AG | [www.bbraun.com/diabetes](http://www.bbraun.com/diabetes)